

PAJUNK®


*Pediatric Anaesthesia
Cannula Systems for Caudal,
Epidural, and Plexus Anaesthesia*



M E D I C A L T E C H N O L O G Y



Pediatric Anaesthesia




Fine and precise:

The PAJUNK® Special Program for Pediatric Anaesthesia

As manufactured of various specialty cannula, PAJUNK® offers a complete product line designed specifically for pediatric anaesthesia in

regional anaesthesia. We offer high quality cannula for caudal, spinal, epidural, and plexus anaesthesia.


Caudal Anaesthesia



SPROTTE® 27 G



The SPROTTE® has initiated a worldwide trend to atraumatic puncture, which has proven itself in pediatric anaesthesia applications as well.




Crawford cannula 23 G, graduated, with depth stop




The Caudal cannula with its cone point geometry according to Crawford has proven to be particularly suitable in pediatric puncture.

Epidural Anaesthesia





PAJUNK® produces two special types of cannulas for Epidural anaesthesia:


With the well proven Tuohy specifically designed by PAJUNK® with a special polish finish




SPROTTE® SPECIAL tip



Both cannula types with Tuohy and SPROTTE® tips are available as EpiLong Pediatric Sets.






Plexus Anaesthesia


UniPlex cannula "Single Shot": (UP/UPS) cannula



The UniPlex cannula were developed by PAJUNK® specifically for the disposable application in the peripheral conductive anaesthesia. These cannula come in many different gages and sizes. PAJUNK® understands the growing trend of this application in pediatric anaesthesia so we offer two versions a SPROTTE® tip and a Facet tip with a shorter lengths and very small gages.




SPROTTE® tip (UPS cannula)




Facet tip (UP cannula)


Plexolong Set according to Meier "Catheter Technique"






The patented Plexus Catheter Set was developed in cooperation with Dr. Meier specifically for the continuous peripheral techniques for regional anaesthesia for children.





Plexolong cannula 21 G



Tuohy tip





Caudal

Product	Size	Item No.	PU
SPROTTE® cannula	22 G x 50 mm	071151-30C	25
SPROTTE® cannula	24 G x 35 mm	001151-30E	25
SPROTTE® cannula	25 G x 35 mm	001151-29E	25
SPROTTE® cannula	27 G x 35 mm	111151-27A	25
Crawford cannula	22 G x 40 mm	1149-12E040	25
Crawford cannula	23 G x 25 mm	1149-12D025	25



Epidural

Product	Size	Catheter	Item No.	PU
Epilong Set I				
Tuohy cannula	19 G x 45 mm	23 G x 45 cm	0131153-50	10
SPROTTE® SPECIAL cannula	19 G x 45 mm	23 G x 45 cm	0131152-50	10
Epidural cannula				
Tuohy cannula	18 G x 50 mm		1150-4K050	25
Tuohy cannula	19 G x 50 mm		1150-4I050	25

X-ray contrast catheter guide
0.2 µm flat filter
Tuohy Borst adapter
6-ml-LOR syringe Luer slip



Plexus

Product	Size	Catheter	Item No.	PU
UniPlex cannula with Facet tip				
UP 3/35	25 G x 35 mm		001153-80	10
UP 3/25	24 G x 25 mm		001153-75	10
UP 3/40	22 G x 40 mm		001153-70	10
UP 3/50	22 G x 50 mm		001153-74	10
UniPlex UPS				
UPS 3/35	25 G x 35 mm		001151-27A	10
UPS 3/50	25 G x 50 mm		001151-29G	10
UPS 3/40	24 G x 40 mm		001151-30G	10
Plexolong-Set with SPROTTE® SPECIAL tip with Tuohy tip				
	21 G x 30 mm	25 G x 50 cm	531151-31B	10
	21 G x 30 mm	25 G x 50 cm	551151-31C	10



NanoLine

Product	Size	Catheter	Item No.	PU
UniPlex-NanoLine Facet Polish				
UP 3/25	24 G x 25 mm		001156-75	10

PAJUNK GmbH
Medical Technology
Karl-Hall-Strasse 1
D-78187 Geisingen/Germany
Telefon +49 (0) 77 04/92 91-0
Telefax +49 (0) 77 04/92 91-6 00
www.pajunk.com

XS200064F 10/05

